

Safe Transportation of Patients on Oxygen Devices In Airborne Precautions

Oxygen Devices:

❖ **Important:** Patient must wear a surgical mask for transport and in close contact with others.




Non-Rebreather Mask



Oxy-Vent with HME (for Trachs)



❖ **Important:** Oxygen devices not recommended during transport:





Safe Transportation of Patients on Oxygen Devices In Airborne Precautions

Non-Invasive Ventilation (NIV)

- ❖ **Important for single limb circuits:** Replace existing exhalation port adaptor and add viral filter.




Switch to blue elbow on mask if connecting to the SERVO.



See Mist Hut Instructions

Manual Ventilation

- ❖ Utilize the Portex® face mask if readily available to improve seal over the patient's face.
- ❖ T-piece resuscitator for newborns (do not add any filters to device).



Transport Ventilator

- ❖ Place HME and viral filters onto circuit (see picture).
- ❖ **Pause or place on standby so flow from ventilator circuit is stopped prior to disconnection.**

